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粪便DNA提取试剂盒 FastDNA Spin Kit For Feces

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快速高效提取粪便DNA，可裂解任何粪便样品。可直接用于PCR等其下游实验，The FastDNA SPIN Kit for Feces Designed for use with the FastPrep Instruments, host cells as well as bacteria, fungi, viruses, protists and other cells present in fecal samples are easily lysed within 40 seconds. These benchtop devices use a unique, optimized motion that provide an extremely quick, efficient and highly reproducible homogenization that surpasses traditional extraction methods.

Samples are placed into 2.0 mL tubes containing Lysing Matrix E, a mixture of ceramic and silica particles designed to

efficiently lyse all stool organisms. Stool samples are pre-treated by washing in a Pre-Lysis Buffer. Homogenization takes place in the presence of MT Buffer and Sodium Phosphate Buffer. The pre-lysing and lysing reagents were carefully developed to protect and solubilize nucleic acids upon cell lysis, facilitate extraction of genomic DNA, and minimize RNA, humic acid or polyphenol contamination. Following lysis, samples are centrifuged to pellet debris and lysing matrix. DNA is purified from the supernatant with a silica-based GENE CLEAN procedure using SPIN filters. Eluted DNA is ready for digestion, electrophoresis, PCR and any other desired application.

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The Fast DNA SPIN Kit for Feces quickly and efficiently isolates PCR-ready genomic DNA from fresh or frozen human and animal stool samples.

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Product Type

Kits & Assays

Components

Lysing Matrix E 50x 2.0 mL tubes, PLS Solution 25 mL, MT Buffer 8 mL, Sodium Phosphate Buffer 100 mL, PPS 25 mL, Binding Matrix 66 mL, Wash Buffer #1 120 mL, Wash Buffer #2 12 mL, TES 20 mL, BBS Gel Loading Dye 200 μ l, SPIN Filter Tubes 50 each, Catch Tubes 50 each

Instrument Compatibility

FastPrep-24; FastPrep 120; FastPrep-96; SuperFastPrep-1

Sample Type

Soil; Feces; Environmental Water; Wastewater; Sludge

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