# NxGen T4 DNA Ligase RNase

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品牌：Lucigen

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Lucigen 30243-2 30281-2 NxGen T4 DNA Ligase RNase Inhibitor 产品说明书

Reliably join double-stranded nucleic acid chains.可靠地连接双链核酸链。

T4 DNA Ligase catalyzes the formation of a phosphodiester bond between the terminal 5' phosphate and 3' hydroxyl groups of duplex DNA or RNA. The enzyme efficiently joins blunt and cohesive ends and repairs single stranded nicks in duplex DNA, RNA or DNA/RNA hybrids.

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2X Rapid Ligation Buffer is composed of 132 mM Tris-HCl, 20 mM MgCl<sub>2</sub>, 2 mM dithiothreitol, 2 mM ATP, 15% PEG, pH 7.6 @ 25 ° C.

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Figure 1: Purity The purity of the T4 DNA ligase is equal or greater than 95% as judged by SDS-polyacrylamide gel electrophoresis with Coomassie blue staining. Molecular weight markers are in Lane 1. T4 DNA ligase was loaded at 5 units (Lane 2) and 10 units (Lane 3). 图1：纯度通过用考马斯蓝染色的SDS-聚丙烯酰胺凝胶电泳判断，T4 DNA连接酶的纯度等于或大于95%。分子量标记在泳道1中。T4 DNA连接酶以5个单位（泳道2）和10个单位（泳道3）加载。

Figure 2: Exonuclease and Endonuclease Activity Assay The T4 DNA ligase is free of detectable of exo- and endonuclease activities, as judged by agarose gel electrophoresis following incubation of 10 units of enzyme for 16 hours at 37 ° C with 1 μg of HindIII-digested DNA (Lanes 3 and 4) and supercoiled pUC19 DNA (Lanes 5 and 6).

6). Molecular weight markers are in Lanes 1 and 2.

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Unit Definition: One Weiss Unit is defined as the amount of enzyme required to convert 1 nmol of <sup>32</sup>P-labeled inorganic pyrophosphate into Norit adsorbable material in 20 minutes at 37 °C, using specified reaction conditions.

Note: 1 Weiss Unit is approximately 67 cohesive end units.

Source: A recombinant E. coli strain carrying the cloned T4 DNA Ligase gene.

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