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YunTek Super VPI is suitable for front-end gasification and sample transfer with pressure liquid samples entering the analysis instrument for analysis. Its design and performance have certain advantages compared with the same kind of

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This product is suitable for sample processing prior to the following analysis：

1. 各类工业用烯烃中微量一氧化碳、二氧化碳的分析，配合气相色谱仪使用。 Analysis of trace amounts of carbon monoxide and carbon dioxide in various kinds of industrial olefins, used in conjunction with gas chromatography.

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Analysis of the composition of liquid hydrocarbons, used in conjunction with gas chromatography.

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Other component analyses, such as sulphide analysis in hydrocarbons.

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This instrument is more convenient than the general liquid sample valve and the analysis results are more accurate.

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All liquid samples are vaporized and transmitted in a constant pressure environment, which overcomes the errors caused by the change of density in the liquid sample during the gasification process, and ensures that the analysis spectrum peaks well.

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The whole transmission pipeline has adjustable thermal function to ensure that the sample does not condense.

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Sampling of liquid and gas samples can be achieved

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The system has its own cleaning function to reduce cross-contamination

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The whole process of inert pipeline and inert treatment filter, to prevent corrosion, adsorption.

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Trace substances that can be used for the detection of sulfides and oxides.

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Has the pressure accurate control function, makes each sample pressure consistent, guarantees the analysis effect.

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Two user-adjustable heaters are designed to control the temperature of the flash steam and the feeding pipe

respectively, so as to achieve the full range of heat and prevent the second condensation after the sample is vaporized.

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After the flash steam sample, repeated sample analysis, the analysis results relative to the standard deviation of the

maximum control at 0.5 %, superior to the same domestic products.

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高气化温度 Maximum gasification temperature

200

温控精度 Temperature control accuracy

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Petrochemical customers 'sample components are dominated by liquid carbon. Theoretically, there is a problem of incomplete volatilisation in the gaseous sample of hot water bath in liquid samples, and the recurrence of volatile components is poor (the repeatability of heavy hydrocarbon components is poor, even greater than 50 %). The actual data also verified that the customer's direct gaseous sample data is poorly repeatable, that is, the use of domestic brand

flash evaporators. Due to the use of excessive transmission lines, the sample can not be completely gasified, which in turn will lead to incomplete sample analysis results and other errors., There is also the phenomenon of insufficient volatilisation or secondary condensation of heavy hydrocarbon components. The peak area of carbonyl sulfur component is less than 0.5 % of the standard difference in the spectrum of sulphide analysis by Flash steaming sample, which ensures the repeatability of the analysis results. The results of YunTek flash evaporator on sulphide composition were proved.

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A gas sample valve on the analyzer interface.

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The analytical instrument needs an external connection port.

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