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Respiratory protective devices, half-masks and surgical masks – An Overview呼吸防护装备，半面口罩和外科手术口罩
interest in the last few months due to the COVID-19 pandemic. Different kinds of masks are increasingly used by people who
available and each is designed for a specific purpose. In the following a short overview of the different protective devices is pr
几个月的热门。人们在公共场合使用越来越多各种不同类型的口罩。人们可以买到各种不同的产品，而这些产
防护设备。We will consider half-masks only which cover the mouth and nose area. Here, basically three different types of r
identified:我们将讨论只遮盖口鼻的半面口罩。在这里，我们可以将口罩大致划分为三个类型：Surgical masks 外
过滤材料制成的半面口罩 Half-mask with exchangeable cassette filters 带可更换盒式过滤器的半面口罩The last two pr
wearer from hazardous substances. Surgical masks on the other hand serve a different purpose. They are designed to protect t

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Lorenz Megertebau's testing machines, such as the FMP03 or FMP04, can measure the filtration efficiency for all particle filters as green in the Table. Additionally, the exhalation and inhalation pressures can be determined as well. They are not stated in the Table, e.g. with new and loaded filters, different flow rates, masks with or without valve, etc. Lorenz Megertebau公司出品的测试机可以测量颗粒过滤器（如石蜡油、塑化剂等）颗粒过滤器的过滤效率，即上表中绿色一栏。此外，还可以确定呼出和吸入压力。上表没有列出呼出和吸入压力，例如，过滤器是否全新、流速、带不带阀门等。More details about e.g. particle size distribution and testing methods are listed below. 有关颗粒大小分布和测试方法的更多细节，请参阅以下内容。Particle filters颗粒过滤器So far particle filtration efficiency is always the same: to remove particles from inhaled air that are potentially capable of entering the human respiratory system. 一直以来，人们一直在谈论颗粒过滤，但是并没有具体规范。颗粒过滤的目的只有一个：去除吸入空气中的颗粒。颗粒大于15 μm 的颗粒会被人体自动过滤。非常小的颗粒，如病毒，直径小于2.5 μm 的颗粒甚至可以进入肺泡。这两种颗粒物分别被称为PM10 (ca. 2.5 – 10 μm diameter) and PM2.5 (< 2.5 μm diameter) respectively. For comparison, bacteria are much larger with diameters between 0.5 μm to 3 μm, while viruses are much smaller with the majority between 0.02 to 0.3 μm. 大的颗粒会被人体自动过滤。非常小的颗粒，如病毒，直径小于2.5 μm 的颗粒甚至可以进入肺泡。这两种颗粒物分别被称为PM10 (ca. 2.5 – 10 μm diameter) and PM2.5 (< 2.5 μm diameter) respectively. For comparison, bacteria are much larger with diameters between 0.5 μm to 3 μm, while viruses are much smaller with the majority between 0.02 to 0.3 μm. Particle filters are tested with a certain concentration and size distribution. Different standards require different substances and slightly different size distributions. Common aerosol generators produce aerosols whose particle diameters are (logarithmically) normally distributed, the well-known Gaussian distribution (normal distribution). The most intuitive way is to state how many particles are present at each diameter (count distribution). The curve's maximum is called the "count median diameter" (CMD). 普通气溶胶发生器生成气溶胶的颗粒大小分布符合对数正态分布，即著名的高斯钟形曲线。其特征是最大位置和宽度（标准偏差）。最直观的方法是描述各直径（或者直径区间）的颗粒数量。最大值被称为“计数中位直径”（CMD）。But there are also other definitions, such as the "mass distribution". When there were no optical particle counters available that we have today. The only option was to determine the amount of particles in a certain distribution was, and still is, of high interest. It is very important not to confuse these two types of distributions and always clearly distinguish between them. Analogously called the "mass median diameter" (MMD). For the same aerosol the MMD is always larger than the CMD. 是，也可以采用其它方式来定义，如“质量分布”。在科学家按照颗粒大小对气溶胶进行分类时，还没有光学颗粒计数器。因此，质量分布仍然广泛使用。请注意，不要将两种分布类型混淆，同时清楚地表明要求的分布类型。在同一种气溶胶来说，质量中位直径总是大于计数中位直径。Testing methods测试方法Filtration efficiencies are measured with a test specimen. Photometers before and after the filter determine the respective aerosol concentrations which yield the efficiency. PAO or DOP. Solid aerosols in the described standards always consist of sodium chloride particles. 令特定气溶胶通过测试确定前后的气溶胶浓度，进而确定过滤效率。液态气溶胶一般包含石蜡油、癸二酸二异辛酯、烯烃或者塑化剂。Inhalation/exhalation resistance is determined with a differential pressure sensor connected to the filter mount. Some devices can determine inhalation/exhalation pressures automatically with integrated pneumatic valves to switch the flow direction. 吸入/呼出阻力由连接到过滤器的压差传感器确定。某些设备可以通过集成的气动阀更改气流方向，实现自动测量吸入和呼出压力。The bacterial filtration efficiency (BFE) is tested with a mean diameter in the range of (3.0 ± 0.3) μm. These droplets contain live Staphylococcus aureus bacteria. Droplets are incubated. The number of colonies compared to that for control runs give the efficiency. No particle filtration test is performed for surgical masks. The bacterial filtration efficiency (BFE) is tested with a mean diameter in the range of (3.0 ± 0.3) μm. These droplets contain live Staphylococcus aureus bacteria. Droplets are incubated. The number of colonies compared to that for control runs give the efficiency. No particle filtration test is performed for surgical masks. 手术口罩的细菌过滤效率。这些飞沫含活性金黄色酿脓葡萄球菌。通过口罩的飞沫将聚集在琼脂板上，同时培养细菌。无需颗粒过滤测试。