

# BIO-RAD伯乐T100 PCR基因扩增仪

|      |                                    |
|------|------------------------------------|
| 产品名称 | BIO-RAD伯乐T100 PCR基因扩增仪             |
| 公司名称 | 北京赛百奥科技有限公司                        |
| 价格   | .00/个                              |
| 规格参数 | 品牌:BIO-RAD<br>型号:186-1096<br>产地:美国 |
| 公司地址 | 北京市海淀区亮甲店130号华玉大厦505               |
| 联系电话 | 010-51282213 13601128645           |

## 产品详情

BIO-RAD 公司具有20多年的PCR仪生产经验，引领业界半导体PCR仪的潮流，满足广大客户的不同需求。T100 PCR仪传承伯乐性能稳定、技术创新的理念，是科研人员的明智之选。流线型简洁设计的T100包括诸多优点，如智能化触屏操作，具有8道编程梯度温度以及快速、稳定的控温性能，让您的PCR实验如此简单。

节省编程时间：直观的触屏操作，使编程更为便捷

快速优化实验流程：温度梯度功能确保在一次实验中优化PCR条件，快速获得佳结果

节省实验室空间：机身小巧紧凑，更坚固，便于安防

先进的文件管理模式：个性化文件夹及USB闪存，易于保存和备份程序

数据的高度一致性：性能稳定，升降温速快，重复性高，结果可靠

## 技术参数：

|                 |   |
|-----------------|---|
| Sample capacity | 96 x 0.2 ml tubes, 0.2 ml tube strips, or 1 x 96-well plate |
|-----------------|---|

|                                |   |
|--------------------------------|---|
| Maximum ramp rate, ° C/sec     | 4   |
| Average ramp rate, ° C/sec     | 2.5   |
| Temperature range              | 4 – 100 ° C   |
| Temperature accuracy           | ± 0.5 ° C of programmed target  |
| Temperature uniformity         | ± 0.5 ° C well-to-well within 30 sec of arrival at target temperature |
| Input power                    | 100 – 150 VAC, 50 – 60 Hz; 220 – 240 VAC, 50 – 60 Hz; 670 W maximum   |
| Display                        | 5.7" VGA color touch screen   |
| Port                           | 1 USB A   |
| Memory                         | 500 typical programs; unlimited with USB flash drive expansion        |
| Dimensions (W x D x H)         | 26 x 47 x 23 cm (10 x 18 x 9")  |
| Weight                         | 9 kg (20 lb)  |
| Gradient accuracy              | ± 0.5 ° C of programmed temperature                                   |
| Row uniformity                 | ± 0.5 ° C well-to-well (within row) within 30 sec                     |
| Gradient range                 | 30 – 100 ° C  |
| Temperature differential range | 1 – 25 ° C  |

订货号：1861096

更多详细资料请登录[www.cbio21.com](http://www.cbio21.com)查询