

## 栓塞弹簧圈 库克cook

产品名称	栓塞弹簧圈 库克cook
公司名称	上海美伯斯医疗科技有限公司
价格	500.00/个
规格参数	品牌:库克cook 型号:imwce 产地:美国
公司地址	中国（上海）自由贸易试验区临港新片区老芦公路536号
联系电话	13770699905 13770699905

## 产品详情

产品名称：库克栓塞弹簧圈、库克弹簧圈、栓塞弹簧圈库克栓塞弹簧圈生产商：库克（中国）医疗贸易有限公司

### 普通弹簧圈0.35系统

MWCE-25-1.2-2, MWCE-25-2.5-2, MWCE-25-2-3, MWCE-25-2-5, MWCE-25-3-2, MWCE-25-3-3, MWCE-25-3-4, MWCE-25-3-5, MWCE-25-4-3, MWCE-25-5-5, MWCE-35-1.2-2, MWCE-35-1.5-2, MWCE-35-10-10, MWCE-35-10-12, MWCE-35-10-15, MWCE-35-10-4, MWCE-35-10-5, MWCE-35-10-8, MWCE-35-1-2, MWCE-35-12-10, MWCE-35-12-5, MWCE-35-1-3, MWCE-35-14-12, MWCE-35-14-5, MWCE-35-14-8, MWCE-35-15-10, MWCE-35-15-12, MWCE-35-15-15, MWCE-35-15-20, MWCE-35-2.5-2, MWCE-35-2-2, MWCE-35-2-3, MWCE-35-2-4, MWCE-35-2-5, MWCE-35-2-8, MWCE-35-3-10, MWCE-35-3-2, MWCE-35-3-3, MWCE-35-3-4, MWCE-35-3-5, MWCE-35-3-6, MWCE-35-3-7, MWCE-35-3-8, MWCE-35-4-2, MWCE-35-4-3, MWCE-35-4-4, MWCE-35-4-6, MWCE-35-4-7, MWCE-35-4-8, MWCE-35-5-10, MWCE-35-5-12, MWCE-35-5-15, MWCE-35-5-2, MWCE-35-5-3, MWCE-35-5-4, MWCE-35-5-5, MWCE-35-5-6, MWCE-35-5-7, MWCE-35-5-8, MWCE-35-6-10, MWCE-35-6-3, MWCE-35-6-4, MWCE-35-6-5, MWCE-35-6-6, MWCE-35-7-3, MWCE-35-7-5, MWCE-35-8-10, MWCE-35-8-12, MWCE-35-8-15, MWCE-35-8-4, MWCE-35-8-5, MWCE-35-8-6, MWCE-35-8-8, MWCE-38-1.2-2, MWCE-38-10-10, MWCE-38-10-12, MWCE-38-10-15, MWCE-38-10-20, MWCE-38-10-5, MWCE-38-10-6, MWCE-38-10-7, MWCE-38-10-8, MWCE-38-12-10, MWCE-38-12-14, MWCE-38-12-5, MWCE-38-12-8, MWCE-38-1-2, MWCE-38-1-3, MWCE-38-13-8, MWCE-38-15-10, MWCE-38-15-12, MWCE-38-15-15, MWCE-38-15-2, MWCE-38-15-20, MWCE-38-15-25, MWCE-38-15-30, MWCE-38-15-8, MWCE-38-16-45, MWCE-38-2.5-2, MWCE-38-2-2, MWCE-38-2-3, MWCE-38-2-5, MWCE-38-2-8, MWCE-38-3-2, MWCE-38-3-3, MWCE-38-3-4, MWCE-38-3-5, MWCE-38-3-8, MWCE-38-4-12, MWCE-38-4-3, MWCE-38-4-4, MWCE-38-4-5, MWCE-38-4-6, MWCE-38-4-7, MWCE-38-5-10, MWCE-38-5-12, MWCE-38-5-15, MWCE-38-5-2, MWCE-38-5-20, MWCE-38-5-3, MWCE-38-5-4, MWCE-38-5-5,

MWCE-38-5-6, MWCE-38-5-7, MWCE-38-5-8, MWCE-38-6.5-5, MWCE-38-6-3, MWCE-38-6-4,  
MWCE-38-6-5, MWCE-38-6-6, MWCE-38-6-8, MWCE-38-7-5, MWCE-38-8-10, MWCE-38-8-12,  
MWCE-38-8-15, MWCE-38-8-3, MWCE-38-8-4, MWCE-38-8-5, MWCE-38-8-6,  
MWCE-38-8-8, MWCE-38-9-5, MWCE-52-10-10, MWCE-52-10-15, MWCE-52-10-6, MWCE-52-10-8,  
MWCE-52-15-10, MWCE-52-15-12, MWCE-52-15-15, MWCE-52-15-20, MWCE-52-5-5, MWCE-52-5-6,  
MWCE-52-5-8, MWCE-52-6-5, MWCE-52-6-6, MWCE-52-8-5, MWCE-52-8-6, MWCE-52-8-8

普通弹簧圈imwce

IMWCE-35-1-2;IMWCE-35-1-3;IMWCE-35-1-5;

IMWCE-35-1.5-3;IMWCE-35-1.5-5;IMWCE-35-2-2;

IMWCE-35-2-3;IMWCE-35-2-4;IMWCE-35-2-5;IMWCE-35-2-8;

IMWCE-35-2.5-5;IMWCE-35-3-2;IMWCE-35-3-3;

IMWCE-35-3-4;IMWCE-35-3-5;IMWCE-35-3-6;IMWCE-35-3-7;

IMWCE-35-3-8;IMWCE-35-3-10;IMWCE-35-4-2;

IMWCE-35-4-3;IMWCE-35-4-4;IMWCE-35-4-5;IMWCE-35-4-6;

IMWCE-35-4-7;IMWCE-35-4-8;IMWCE-35-5-2;IMWCE-35-5-3;

IMWCE-35-5-4;IMWCE-35-5-5;IMWCE-35-5-6;IMWCE-35-5-7;

IMWCE-35-5-8;IMWCE-35-5-10;IMWCE-35-5-12;

IMWCE-35-5-15;IMWCE-35-6-3;IMWCE-35-6-4;

IMWCE-35-6-5;IMWCE-35-6-6;IMWCE-35-6-7;IMWCE-35-6-8;

IMWCE-35-6-10;IMWCE-35-7-3;IMWCE-35-7-5;

IMWCE-35-8-4;IMWCE-35-8-5;IMWCE-35-8-6;IMWCE-35-8-8;

IMWCE-35-8-10;IMWCE-35-8-12;IMWCE-35-8-15;

IMWCE-35-10-4;IMWCE-35-10-5;IMWCE-35-10-6;

IMWCE-35-10-8;IMWCE-35-10-10;IMWCE-35-10-12;

IMWCE-35-10-15;IMWCE-35-10-20;IMWCE-35-12-5;

IMWCE-35-12-8;IMWCE-35-12-10;IMWCE-35-14-3;

IMWCE-35-14-5;IMWCE-35-14-6;IMWCE-35-14-8;

IMWCE-35-14-12;IMWCE-35-15-5;IMWCE-35-15-10;

IMWCE-35-15-12;IMWCE-35-15-15;IMWCE-35-15-20;

IMWCE-35-20-10;IMWCE-35-20-15;IMWCE-35-20-20;

IMWCE-38-1-2;IMWCE-38-5-2;IMWCE-38-5-3;IMWCE-38-5-4;

IMWCE-38-5-5;IMWCE-38-5-6;IMWCE-38-5-7;IMWCE-38-5-8;

IMWCE-38-5-9;IMWCE-38-5-10;IMWCE-38-5-12;

IMWCE-38-5-15;IMWCE-38-5-20;IMWCE-38-6-10;

IMWCE-38-6-12;IMWCE-38-6-15;IMWCE-38-8-3;

IMWCE-38-8-4;IMWCE-38-8-5;IMWCE-38-8-6;IMWCE-38-8-8;

IMWCE-38-8-10;IMWCE-38-8-12;IMWCE-38-8-15;

IMWCE-38-10-8;IMWCE-38-10-10;IMWCE-38-10-12;

IMWCE-38-10-15;IMWCE-38-10-20;IMWCE-38-12-14;

IMWCE-38-15-10;IMWCE-38-15-12;IMWCE-38-15-15;

IMWCE-38-15-20;IMWCE-38-15-25;IMWCE-38-15-30;

IMWCE-38-16-45;IMWCE-38-20-10;IMWCE-38-20-15;

IMWCE-38-20-20;

#### TORNADO塔型弹簧圈

MWCE-18S-3/2-TORNADO;MWCE-18S-4/2-TORNADO;

MWCE-18S-4/3-TORNADO;MWCE-18S-5/2-TORNADO;

MWCE-18S-6/2-TORNADO;MWCE-18S-7/2-TORNADO;

MWCE-18S-7/3-TORNADO;MWCE-18S-8/4-TORNADO;

MWCE-18S-8/5-TORNADO;MWCE-18S-10/4-TORNADO;

MWCE-35-10/3-TORNADO;MWCE-35-10/4-TORNADO;

MWCE-35-10/5-TORNADO;MWCE-35-4/3-TORNADO;

MWCE-35-5/3-TORNADO;MWCE-35-6/3-TORNADO;

MWCE-35-7/3-TORNADO;MWCE-35-8/4-TORNADO;

MWCE-35-8/5-TORNADO;MWCE-38-4/3-TORNADO;

MWCE-38-7/3-TORNADO;MWCE-38-10/4-TORNADO;

MWCE-38-10/5-TORNADO;MWCE-18S-4/8-TORNADO-LEF;

MWCE-18-5-5-NESTER;MWCE-18-7-5-NESTER;

MWCE-18-14-5-NESTER;MWCE-18-14-4-NESTER;

MWCE-18-14-6-NESTER;MWCE-18-14-8-NESTER;

MWCE-18-14-10-NESTER;MWCE-35-7-3-NESTER;

MWCE-35-7-5-NESTER;MWCE-35-14-3-NESTER;

MWCE-35-14-5-NESTER;MWCE-35-14-4-NESTER;

MWCE-35-14-6-NESTER;MWCE-35-14-8-NESTER;

MWCE-35-14-10-NESTER;MWCE-35-14-12-NESTER;

MWCE-35-20-5-NESTER;MWCE-35-20-6-NESTER;

MWCE-35-20-8-NESTER;MWCE-35-20-10-NESTER;

MWCE-35-20-12-NESTER;MWCE-35-20-14-NESTER;

MWCE-35-20-16-NESTER;MWCE-38-7-3-NESTER;

MWCE-38-7-5-NESTER;MWCE-38-14-3-NESTER;

MWCE-38-14-5-NESTER;MWCE-38-20-5-NESTER;

MWCE-38-20-6-NESTER;MWCE-38-20-8-NESTER;

MWCE-38-20-10-NESTER;MWCE-38-20-12-NESTER;

MWCE-38-20-14-NESTER;MWCE-38-20-18-NESTER;

HILAL直筒微弹簧圈

MWCE-18S-1.2-2-A

MWCE-18S-1.8-3-A

MWCE-18S-2.4-4-A

MWCE-18S-2.0-5-B

MWCE-18S-3.0-7-B

MWCE-18S-4.2-10-B

MWCE-18S-5.2-12-B

MWCE-18S-2.6-3-C

MWCE-18S-3.5-4-C

MWCE-18S-0.5-0-HILAL

MWCE-18S-0.7-0-HILAL

MWCE-18S-1.0-0-HILAL

MWCE-18S-1.0-3-HILAL

MWCE-18S-1.5-0-HILAL

MWCE-18S-1.5-5-HILAL

MWCE-18S-2.0-2-HILAL

MWCE-18S-2.0-4-HILAL

MWCE-18S-2.1-7-HILAL

MWCE-18S-3.0-3-HILAL

MWCE-18S-3.0-4-HILAL

MWCE-18S-3.0-10-HILAL

MWCE-18S-4.0-6-HILAL

MWCE-18S-4.0-7-HILAL

MWCE-18S-6.0-5-HILAL

MWCE-18S-6.0-7-HILAL

MREYE核磁圈 (MRI款)

NESTER鸟巢圈

MWCE-18-14-10-NESTER, MWCE-18-14-3-NESTER, MWCE-18-14-4-NESTER, MWCE-18-14-5-NESTER,  
MWCE-18-14-6-NESTER, MWCE-18-14-8-NESTER, MWCE-18-3-2-NESTER, MWCE-18-3-3-NESTER,  
MWCE-18-4-3-NESTER, MWCE-18-5-2-NESTER, MWCE-18-5-3-NESTER, MWCE-18-5-5-NESTER,  
MWCE-18-7-10-NESTER, MWCE-18-7-2-NESTER, MWCE-18-7-3-NESTER, MWCE-18-7-4-NESTER,  
MWCE-18-7-5-NESTER, MWCE-18-7-6-NESTER, MWCE-18-7-8-NESTER,  
MWCE-35-14-10-NESTER, MWCE-35-14-12-NESTER, MWCE-35-14-14-NESTER, MWCE-35-14-16-NESTER,  
MWCE-35-14-18-NESTER, MWCE-35-14-20-NESTER, MWCE-35-14-3-NESTER, MWCE-35-14-4-NESTER,  
MWCE-35-14-5-NESTER, MWCE-35-14-6, MWCE-35-14-6-NESTER, MWCE-35-14-8-NESTER,  
MWCE-35-20-10-NESTER, MWCE-35-20-12-NESTER, MWCE-35-20-14-NESTER,

MWCE-35-20-16-NESTER,MWCE-35-20-18-NESTER, MWCE-35-20-20-NESTER, MWCE-35-20-5-NESTER, MWCE-35-20-6-NESTER, MWCE-35-20-8-NESTER,MWCE-35-4-4-NESTER, MWCE-35-4-5, MWCE-35-4-5-NESTER,MWCE-35-6-6-NESTER, MWCE-35-7-10-NESTER, MWCE-35-7-12-NESTER,MWCE-35-7-14-NESTER, MWCE-35-7-16-NESTER, MWCE-35-7-18-NESTER, MWCE-35-7-20-NESTER, MWCE-35-7-3-NESTER,MWCE-35-7-4-NESTER, MWCE-35-7-5-NESTER, MWCE-35-7-6-NESTER, MWCE-35-7-8-NESTER, MWCE-35-8-8-NESTER,MWCE-38-14-10-NESTER, MWCE-38-14-12-NESTER, MWCE-38-14-4-NESTER, MWCE-38-14-6-NESTER, MWCE-38-14-8-NESTER,MWCE-38-4-4-NESTER, MWCE-38-4-6-NESTER, MWCE-38-6-6-NESTER, MWCE-38-8-8-NESTER

产品型号：MWCE-25-1.2-2 MWCE-35-1.2-2 MWCE-35-4-2 MWCE-38-1.2-2 MWCE-38-4-12  
MWCE-25-2.5-2 MWCE-35-1.5-2 MWCE-35-4-3 MWCE-38-10-10 MWCE-38-4-3 MWCE-25-2-3  
MWCE-35-10-10 MWCE-35-4-4 MWCE-38-10-12 MWCE-38-4-4 MWCE-25-2-5 MWCE-35-10-12  
MWCE-35-4-5 MWCE-38-10-15 MWCE-38-4-5 MWCE-25-3-2 MWCE-35-10-15 MWCE-35-4-6  
MWCE-38-10-20 MWCE-38-4-6 MWCE-25-3-3 MWCE-35-10-4 MWCE-35-4-7 MWCE-38-10-5  
MWCE-38-4-7 MWCE-25-3-4 MWCE-35-10-5 MWCE-35-4-8 MWCE-38-10-6 MWCE-38-5-10 MWCE-25-3-5  
MWCE-35-10-8 MWCE-35-5-10 MWCE-38-10-7 MWCE-38-5-12 MWCE-25-4-3 MWCE-35-1-2  
MWCE-35-5-12 MWCE-38-10-8 MWCE-38-5-15 MWCE-25-5-5 MWCE-35-12-10 MWCE-35-5-15  
MWCE-38-12-10 MWCE-38-5-2 MWCE-52-10-10 MWCE-35-12-5 MWCE-35-5-2 MWCE-38-12-14  
MWCE-38-5-20 MWCE-52-10-15 MWCE-35-1-3 MWCE-35-5-3 MWCE-38-12-5 MWCE-38-5-3  
MWCE-52-10-6 MWCE-35-14-12 MWCE-35-5-4 MWCE-38-12-8 MWCE-38-5-4 MWCE-52-10-8  
MWCE-35-14-5 MWCE-35-5-5 MWCE-38-1-3 MWCE-38-5-5 MWCE-52-15-10 MWCE-35-14-6  
MWCE-35-5-6 MWCE-38-13-8 MWCE-38-5-6 MWCE-52-15-12 MWCE-35-14-8 MWCE-35-5-7  
MWCE-38-15-10 MWCE-38-5-7 MWCE-52-15-15 MWCE-35-15-10 MWCE-35-5-8 MWCE-38-15-12  
MWCE-38-5-8 MWCE-52-15-20 MWCE-35-15-12 MWCE-35-6-10 MWCE-38-15-15 MWCE-38-6.5-5  
MWCE-52-5-5 MWCE-35-15-15 MWCE-35-6-3 MWCE-38-15-2 MWCE-38-6-3 MWCE-52-5-6  
MWCE-35-15-20 MWCE-35-6-4 MWCE-38-15-20 MWCE-38-6-4 MWCE-52-5-8 MWCE-35-2.5-2  
MWCE-35-6-5 MWCE-38-15-25 MWCE-38-6-5 MWCE-52-6-5 MWCE-35-2-2 MWCE-35-6-6  
MWCE-38-15-30 MWCE-38-6-6 MWCE-52-6-6 MWCE-35-2-3 MWCE-35-7-3 MWCE-38-15-8 MWCE-38-6-8  
MWCE-52-8-5 MWCE-35-2-4 MWCE-35-7-5 MWCE-38-16-45 MWCE-38-7-5 MWCE-52-8-6 MWCE-35-2-5  
MWCE-35-8-10 MWCE-38-2.5-2 MWCE-38-8-10 MWCE-52-8-8 MWCE-35-2-8 MWCE-35-8-12  
MWCE-38-2-2 MWCE-38-8-12 MWCE-35-3-6 MWCE-35-3-10 MWCE-35-8-15 MWCE-38-2-3  
MWCE-38-8-15 MWCE-35-3-7 MWCE-35-3-2 MWCE-35-8-4 MWCE-38-2-5 MWCE-38-8-3 MWCE-35-3-8  
MWCE-35-3-3 MWCE-35-8-5 MWCE-38-2-8 MWCE-38-8-4 MWCE-38-3-5 MWCE-35-3-4 MWCE-35-8-6  
MWCE-38-3-2 MWCE-38-8-5 MWCE-38-3-8 MWCE-35-3-5 MWCE-35-8-8 MWCE-38-3-3 MWCE-38-8-6  
MWCE-38-9-5 MWCE-38-3-4 MWCE-38-8

医疗器械注册证编号：国械注进20173771849

适用范围：该产品被设计用于动静脉血管栓塞手术。

结构及组成：该产品由弹簧圈和装载筒组成。弹簧圈由304不锈钢丝和人造纤维(尼龙66)构成。装载筒材料为304不锈钢。环氧乙烷灭菌，产品一次性使用。货架有效期五年。