

X光洗片机HR-430R医用洗片机

产品名称	X光洗片机HR-430R医用洗片机
公司名称	泰兴市红日医疗设备厂
价格	.00/个
规格参数	
公司地址	泰兴市通江路海关北街28号
联系电话	86 0523 87638844 13901438084

产品详情

最大冲洗宽度17英寸 最小冲洗宽度3x4英寸 冲洗胶片品种x射线、ct、mri等医用胶片 冲洗程序显影、定影(无级) 冲洗量/小时(显影20秒时)10x12 :240张14x17 180张(理论冲量) 温度调节药液:20~40 可调烘干:40~60 药液补充方法自动和手动。每累计冲洗0.5m²补充一次

水洗控制由电磁阀自动控制,无片冲洗自动停止 药液补充调节40~400ml可调 水管接口直径 15mm水量2~2.8kw 体积940x650x1160(长x宽x高)

智能化电脑控制系统	数字显示键盘人机对话方式,可选择标准和非标
深槽式滚轴输片系统	冲洗速度可自由调节控制,并具有超温报警和自
可靠的进片检测系统	装药量大,费用低,送片冲洗平衡、可靠
故障自检报警系统	本机采用我厂自行研制的进片检测系统,其进片
节能和自动停机系统	机器若出现故障能自动显示、报警。
水洗自动控制系统	不连续进片自动切断主电板的节能装置:软片冲
抗氧化、防结晶系统	同时马达自动停止转动,机器进入节能状态,烘
药液过滤系统	机器工作即自动进水,停机则不进水,有效节约
自动补液系统	在预备和待机状态时机器自动清洗滚轴,防止结
	药液自动循环过滤,保证冲洗质量
	自动进片补充及液面补充同时独立操作。
intelligent computer controlled system	digital display, keyset man-machine conversation, ca
	developing method, every solutions temp can be used
deep tank roller film entering system	developing speed can be controlled at option. when th
	give the alarm and shut the power automatically.
credible entering films checking system	install plentiful solution, low fees, films entering system
	the machine is used the advanced films entering check
give an alarm automatically system	ourselves, the films entering checking is very credible
saving and stopping auto system	when the machine meets the troubles, it will display an
water developing controlled auto system	stop films the machine will cut the saving system auto
	heated will delay 2 min then stop heating, fan will del
anti-ox, anti-crystal system	motor will stop moving auto, the machine will be in t
filtrate solution system	by the advanced heater, which have high efficient, so c

	supply solution auto	when the machine is working ,it will be in the enterin
		working , it won ' t enter the water; it can be saved th
		when the machine is starting , it will clean the rollers a
		solutions are filtrated auto to auarantee the developin
		supply the films auto, supply timely and supply the sc

批发说明