

biorad伯乐双向电泳槽

产品名称	biorad伯乐双向电泳槽
公司名称	北京赛百奥科技有限公司
价格	.00/个
规格参数	品牌:BIO-RAD 型号:2-D 产地:美国
公司地址	北京市海淀区亮甲店130号华玉大厦505
联系电话	010-51282213 13601128645

产品详情

PROTEAN II xi 2-D 双向电泳槽

【介绍】

The PROTEAN II xi 2-D 电泳槽包含聚丙烯酰胺管式双向电泳需要的所有组件。

在PROTEAN II xi 2-D 电泳槽内，可同时进行向的管式凝胶电泳和第二向的板式凝胶电泳

运用管式凝胶适配器，一次电泳可聚焦多达16根向IEF管式凝胶

斜面玻板，易于放置管式凝胶的板式凝胶表面

PROTEAN II xi 2-D 电泳槽组件：

管式凝胶中央制冷芯

带密封垫条的电泳中央制冷芯

缓冲液槽和带有电线的盖子

板式凝胶灌胶架、玻璃板和三明治夹

垫圈和塞子

玻璃管

【订货信息】

PROTEAN II xi 2-D Cell

165-1931	PROTEAN II xi 2-D Cell, includes a central cooling core with gaskets, lower buffer chamber, lid with 16 cm diameter of glass plates (with beveled inner plates), four 16 cm sandwich clamps, twenty-four 180 mm long glass tubes (180 mm diameter = spacer thickness), 2 tube gel adaptors, 16 grommets, 16 stoppers, two 1.0 mm 2-D combs, upper buffer dam, casting stand with gaskets, leveling bubble, instructions
165-1932	PROTEAN II xi 2-D Cell, includes a central cooling core with gaskets, lower buffer chamber, lid with 16 cm diameter of glass plates (with beveled inner plates), four 16 cm sandwich clamps, twenty-four 180 mm long glass tubes (180 mm diameter = spacer thickness), 2 tube gel adaptors, 16 grommets, 16 stoppers, two 1.5 mm 2-D combs, upper buffer dam, casting stand with gaskets, leveling bubble, instructions
165-1933	PROTEAN II xi 2-D Cell, includes a central cooling core with gaskets, lower buffer chamber, lid with 20 cm diameter of glass plates (with beveled inner plates), four 20 cm sandwich clamps, twenty-four 180 mm long glass tubes (180 mm diameter = spacer thickness), 2 tube gel adaptors, 16 grommets, 16 stoppers, two 1.0 mm 2-D combs, upper buffer dam, casting stand with gaskets, leveling bubble, instructions
165-1934	PROTEAN II xi 2-D Cell, includes a central cooling core with gaskets, lower buffer chamber, lid with 20 cm diameter of glass plates (with beveled inner plates), four 20 cm sandwich clamps, twenty-four 180 mm long glass tubes (180 mm diameter = spacer thickness), 2 tube gel adaptors, 16 grommets, 16 stoppers, two 1.5 mm 2-D combs, upper buffer dam, casting stand with gaskets, leveling bubble, instructions
	Note: Each PROTEAN II xi 2-D cell includes a central cooling core with gaskets, lower buffer chamber with cables, 2 sets of glass plates (with beveled inner plates), 4 sandwich clamps, 24 glass tubes (180 mm diameter = spacer thickness), 2 tube gel adaptors, 16 grommets, 16 stoppers, two 2-D combs (2.0 mm diameter), upper buffer dam, casting stand with gaskets, leveling bubble, instructions. Each cell includes sandwich clamps for running gels of the corresponding size (either 16 cm or 20 cm); 1.0 mm and 1.5 mm indicate the spacer and combs included with the cell.

Tube Gel Casting Accessories

165-2020	Model 225 Tube Gel Casting Stand
165-3122	5.0 mm ID Glass Tubes, 7.0 mm OD, 125 mm length, 24
165-3136	1.0 mm ID Glass Tubes, 6.0 mm OD, 180 mm length, 24
165-3137	1.5 mm ID Glass Tubes, 7.5 mm OD, 150 mm length, 24
165-3138	1.5 mm ID Glass Tubes, 7.5 mm OD, 180 mm length, 24
165-3150	3.4 mm ID Glass Tubes, 5.0 mm OD, 125 mm length, 24
165-3155	2.4 mm ID Glass Tubes, 4.0 mm OD, 160 mm length, 24

Slab Gel Handcasting Accessories

165-1821	Inner Plates, for 16 cm cell, 16 x 20 cm, 2, for PROTEAN II xi cell only
165-1822	Outer Plates, for 16 cm cell, 18.3 x 20 cm, 2, for PROTEAN II xi cell only
165-1823	Inner Plates, for 20 cm cell, 20 x 20 cm, 2
165-1824	Outer Plates, for 20 cm cell, 22.3 x 20 cm, 2
165-1827	Beveled Inner Plates, for 2-D tube gel procedures, 16 cm bevel length, 16 x 20 cm, 2, for PROTEAN II xi cell only
165-1828	Beveled Inner Plates, for 2-D tube gel procedures, 16 cm bevel length, 20 x 20 cm, 2, for PROTEAN II xi cell only

165-1843	PROTEAN II xi Spacers, 16 cm, 1.0 mm spacers
165-1848	PROTEAN II xi Spacers, 20 cm, 1.0 mm spacers
165-1897	PROTEAN II xi Comb, 2-D (1 reference well), 1.0 mm
165-1901	PROTEAN II Sandwich Clamps, 16 cm set (1 left, 1 right)
165-1902	PROTEAN II Sandwich Clamps, 20 cm set (1 left, 1 right)
165-1911	Slab Gel Casting Stand, with gaskets
Other Accessories	
165-1806	Central Cooling Core, with gaskets
165-1807	Buffer Tank
165-1909	Upper Buffer Dam
165-1943	Tube Gel Loading Needles, 18 cm, 22 gauge, blunt tip, luer hub (for casting monomer in small-diam)
165-1944	Tube Gel Extrusion Needles, 9 cm, 26 gauge, beveled tip, luer hub (for removing gels from tubes), 2
System Accessories	
165-2025	PROTEAN II xi Multi-Gel Casting Chamber, includes casting chamber, sealing plate, silicone gasket, separation sheets, 4 acrylic blocks, 10 xi alignment cards, tapered luer connector, leveling bubble, ins
165-4121	Model 495 Gradient Former, 100 – 1,500 ml, includes body with valve stem and tubing connection k

更多详细资料请登录www.cbio21.com查询