一次性塑料试管

产品名称	一次性塑料试管
公司名称	江苏康捷医疗器械有限公司
价格	.00/普通
规格参数	
公司地址	泰州江苏姜堰市人民北路刁家桥
联系电话	0523-88660056 13961082327

产品详情

email: jskj@js-kj.com

一次性使用真空采血管主要由采血管管体,丁基胶内塞,纹面塑料安全帽,抗凝剂或凝固促进剂或血清分离剂组成,抗凝剂分为edta-k2,edta-k3,edta-na2,枸橼酸钠,肝素钠,肝素锂或草酸钾/氟化物,edta/氟化物,肝素/氟化物,sps溶液氟化钠,枸橼酸钠葡萄糖,凝固促凝剂为玻璃粉,附加物:血清分离剂为-烯烃与无水马来酸的共聚物。按照无菌状态分为无菌类和非无菌类两种,无菌采血管试管材质为pet,非无菌采血管试管材质为玻璃或pet;按附加剂和采血容量不同分为若干规格。公称采血量:1ml-10ml。无菌采血管应无菌。 jiangsu kangjie medical devices co,ltd. isfounded in 2004.we are a moderm enterprise. we mainly producemedical consumable goods and lab equipments. for several yearsrapid development, kangjie has been a bright star in inspectionindustry. we have advanced equipments. we have a purification workshopwhich conforms to gmp standard. we have automatic injectionworkshop、standardized production line、quality testcenter、etc. we adhere to "live on quality, develop on technology". weestablished an exacting quality control system according to theinternational standard. we have got ce certificate. such measuresassured our products' quality stability and reliability. thesale web of our products has reached the whole nation andoversea. products: vacuum blood collecting tubes、tip tube swabseries petri dish series、cuvette cup series、sample cup series 地址:

江苏省姜堰市人民北路刁家桥电话:0523-88660056、88660201、88660202传真:0523-88661197邮箱:jskj@js-kj.com网址:http://www.js-kj.comjiangsu kangjie medical devices co,ltd. address: shengao town people

village, jiangyan city, jiangsu,china 225500 tel: 86 523 88660056 fax: 86 523 88661197 website: http://www.js-kj.com