

Disposable Trocar/Disposable Sheath Trocar/Trocar

产品名称	Disposable Trocar/Disposable Sheath Trocar/Trocar
公司名称	江苏新智源医学科技有限公司
价格	面议
规格参数	卫生许可证:苏食药监械生产许2010-0090号 灭菌方式:环氧乙烷 有效期:24
公司地址	泰州市中国医药城三期
联系电话	15161063071

产品详情

types : a/b/e

intended use : the product is indicated for use to establish a port of entry for endoscopic instruments in patients requiring minimally invasive surgical procedures.

efficient — integrated general cannula cap

no need exchange adapter during surgery.

save operation time.

convenient - detachable design , easy for specimen retrieval

detachable cannula cap , easy for specimen and tissue retrieval during surgery .

stable -double sealing , prevent gas leakage

special material , durable double sealing system .

reduce gas leakage when instruments get in or out of the trocar .

maintain artificial pneumoperitoneum even instruments get in or out of the trocar for many times .

stable gas tightness meet the demand of high-level laparoscopic surgery .

reliable- greater abdominal wall retention

the cannula's integrated thread design provides greater abdominal wall retention .

minimal trocar slip-outs .

efficient -low system drag force

low system drag force and funnel lead-in .

enables smooth, one-handed instrument exchange.

innovative-bladeless , optical tip

direct visualization of tissue layers eliminates blind entry insertion.

the bladeless tip separates, rather than cuts, along tissue fibers, pushing tissue and vessels away. reduce trauma, facilitate wound healing .

various -multiple specifications , different materials for choice

various specifications meet different demands from clinical surgeries .

stainless steel or medical polymer materials trocars for different choices

本产品的卫生许可证是苏食药监械生产许2010-0090号，灭菌方式是环氧乙烷，有效期是24，一次性使用是是，药（械）准字是苏食药监械（准）字2011第2080857号