

TwistAmp fpg试剂盒

产品名称	TwistAmp fpg试剂盒
公司名称	畅宇云商（北京）科技有限公司
价格	面议
规格参数	品牌:TwistDX 型号:fpg 货号:TAFPG01kit_T
公司地址	北京市门头沟区石龙经济开发区永安路20号3号楼A-7427室
联系电话	13810792952 13810792952

产品详情

TAFPG01kit_T	TwistAmp™ fpg	用于real-time探针检测
--------------	---------------	-----------------

The fpg kit is tailored for users who want to combine the TwistDx amplification technology with an alternative TwistDx reporter probe system, fluorescent TwistAmp® fpg Probe, in a homogenous format.

Our new formulation kit includes extra ssDNA binding protein and dNTP to improve RPA performance. We highly recommend customers switch to this formulation as we expect improved performance for most uses.

Features:

NEW IMPROVED FORMULATION

Real-time probe detection

Single-molecule detection

Twista® - compatible

Multiplexable

1 x box of 96 reactions

PLEASE NOTE:

For magnetic mixing using the T16 device, bearing balls may be added to reactions at any time before the reaction is initiated with Mg.

These probes have fewer design constraints than TwistAmp®exo Probe, although kinetics of fluorescence accumulation may be slower. In addition to the basic components it contains a powerful nuclease (fpg) which will process TwistAmp fpg Probe during the amplification reaction itself and generate a real-time read-out. The presence of fpg will not reduce the final overall yield of amplified material in contrast to the use of Exonuclease III in the TwistAmp®exo kit allowing end point gel analysis, but the kinetics of probe cutting are sometimes slower than with the exo system.

There are 96 reactions in each kit.