

厂家现货热销流感检测用受体破坏酶RDE 340016 RDE(II) 340122

产品名称	厂家现货热销流感检测用受体破坏酶RDE 340016 RDE(II) 340122
公司名称	武汉御冷科技有限公司
价格	.00/个
规格参数	
公司地址	武汉市东西湖区新沟镇街道办事处卧龙路81号变频节能风机生产项目4号厂房2层（注册地址）
联系电话	17396110611

产品详情

Denka Seiken 日本生研受体破坏酶 340016 (370013) 340122

Denka Seiken 340016 (370013) RDE (受体破坏酶)

20ml*5

Denka Seiken 340122 RDE (受体破坏酶) 20ml*5

贮藏方法/有効期间 贮藏方法 2～10℃ で保存

有効期间 1年(ラベルに表示の使用期限内にご使用ください。)

包装単位 各10mL(液状) 1本

インフルエンザウイルスHI試薬「生研」A型HA抗原*4

インフルエンザウイルスHI試薬「生研」B型HA抗原*4

当社別売として、次の製品も用意してありますのでご利用ください。

各1mL分(凍結乾燥)1本

インフルエンザウイルスHI試薬「生研」A型HI抗血清*4

インフルエンザウイルスHI試薬「生研」B型HI抗血清*4

*4：商品番号は、インフルエンザウイルスの株によって変わりますので、別途お問い合わせください。

各20mL分(凍結乾燥)5本

RDE「生研」(商品番号340016)

RDE(II)「生研」(商品番号340122)

Supplemental Reagents for HI Test

Cat No.	Description	Specification
340-016	RDE	20ml x 5
340-122	RDE (II)	20ml x 5
340-023	DGV for Rubella Virus HI Test	50ml x 10
340-030	PBS(-) for Rubella Virus HI Test	10ml x 10
340-047	25% Kaolin Suspension for Rubella Virus HI Test	20ml x 10
340-054	Diluent for Rubella Virus HI Test	20ml x 10
320-063	14% Kaolin Suspension for R-HI	40ml x 10
320-032	Diluent for Rubella Virus R-HI	45ml x 10
320-049	2% RBC Suspension for R-HI	4ml x 5
320-056	50% RBC Suspension for R-HI	3ml x 5
340-085		

PBS for Measles Virus HI Test 50ml x 5 340-092 25% Kaolin Suspension for Measles HI
Test 20ml x 5 340-108 Diluent for Measles Virus HI Test 20ml x 5 Supplemental
Reagents fo CF Test Cat No. Description Specification 430-014 Complement(Lyophilized)
1ml x 10 430-021 Complement(Lyophilized) 5ml x 5 430-076 Heomlysin
3ml 430-069 Greenss solution :SEIKEN" 20ml Others Cat No. Description
Specification 340-061 Trypsin for Cell Dispersion 150ml 340-078 Trypsin-EDTA
150ml 430-038 Sheep Blood Cells 50ml