

德立功FY15810-D注射器密合性负压测试仪

产品名称	德立功FY15810-D注射器密合性负压测试仪
公司名称	昆山德立功自动化设备有限公司
价格	.00/个
规格参数	品牌:德立功
公司地址	昆山市花桥镇立德企业家园区1号楼1室（注册地址）
联系电话	18511156696

产品详情

FY15810-D注射器密合性负压测试仪

Negative pressure tester for syringe sealing

采用5.7英寸（320 × 240）多彩液晶触摸显示屏，中文菜单显示。并由键盘进行选择，并可由机载打印机打印出测试结果。完全符合GB15810-2001《一次性使用无菌注射器》标准、ISO7886-1-2016标准中的相关标准设计制造。

负压检测：88Kpa误差不大于 ± 0.5Kpa、

负压输出：88Kpa ± 0.5Kpa;

测试时间：1S ~ 10min任意设定，（测量系统负压保持5分钟内不漂移 ± 0.5Kpa）

净重：8.5kg 外形尺寸：400 × 290 × 240（mm）

部分清单：

1.主机一台 2.打印纸1卷 3.质保一年 4.操作说明

5.仪器终身免费依国家标准规定升级服务

A 5.7 inch (320 x 240) multi-color LCD touch screen is used to

display the Chinese menu. It can be selected by keyboard and

printed out by airborne printer. Fully conform to GB15810-2001

"disposable sterile syringe" standards, ISO7886-1-2016

standards in the design and manufacture of relevant standards.

Negative pressure detection: 88 Kpa error is not more than

(+0.5 Kpa)

Negative pressure output: 88 Kpa + 0.5 Kpa;

Testing time: 1 S-10 min arbitrarily set, (the negative

pressure of the measuring system is maintained for 5 minutes

without drift (+0.5 Kpa)

Net weight: 8.5Kg shape size: 400 *290 *240 (mm)

Part of the list:

1. Host 1. 2. Printing paper 1 roll 3. Quality assurance 1

year 4. Operational instructions

5. Instrument upgrade service free of charge for life in

accordance with national standards